



Canine Lymphoma Blood Test



<p>Serum Collection Protocol</p> <ul style="list-style-type: none"> Take 1ml of venous blood and place in a red top tube Allow to clot for 20 minutes at room temperature Centrifuge the sample to separate serum Decant serum and keep refrigerated until ready to ship Complete the Avacta submission form Wrap serum sample in tissue and seal in a biohazard bag. Place in a padded envelope with a small ice pack, and send to Avacta at the address below using FedEx (or equivalent) standard overnight 	<p>Veterinarian:</p> <p>Practice name and address:</p> <p>Zip:</p> <p>Telephone:</p> <p>Email address (For results)</p>
<p>Animal's name:</p> <p>Owner's name:</p> <p>Age: Sex: Neutered: Y / N</p> <p>Breed: Date sampled:</p>	

REQUIRED INFORMATION FOR ALL TESTS:

This test is for: (please circle) MONITORING or DIAGNOSIS

Lymphadenopathy currently: (please circle) YES/ NO

Please Note: Results may be unrepresentative if steroids, NSAIDS, joint supplements or other anti-inflammatory drugs have been administered recently.

FOR CHEMOTHERAPY MONITORING ONLY:

Start date for initial (1st round) chemotherapy _____

Chemotherapy protocol/ drugs used:

Last treatment date: _____ or any other current treatment:

Have you run a previous cLBT for this client? YES/ NO

Clinical History/ Diagnosis (if known): (please print)

To be completed by the owner or on behalf of the owner where possible

I hereby give my consent for the use of this blood sample in the development of the above test and other such tests to diagnose important canine diseases:

Name Signature

Thank you for your support.

Mail samples to:
 Avacta Animal Health, 4011 Discovery Drive,
 Ste N110, Columbia, Missouri, 65211
 Tel: (310) 292-0135 Email: tariq.shah@avacta.com

*Samples submitted are subject to Avacta Animal Health's Terms and Conditions of Business. A copy can be found at:
www.avactaanimalhealth.com/tandc.